# API 20 NE

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3. *Pseudomonas aeruginosa* ATCC 27853
4. *Alcaligenes faecalis* ATCC 35655

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| TESTS           | ACTIVE INGREDIENTS                             | QTY<br>(mg/cup.) | REACTIONS / ENZYMES  | RESULTS                           |                              |
|-----------------|--|------------------|--|-----------------------------------|------------------------------|
|                 |  |                  |  | NEGATIVE                          | POSITIVE                     |
| NO <sub>3</sub> | potassium nitrate                              | 0.136            | reduction of nitrates to nitrites  | NIT 1 + NIT 2 / 5 min             |                              |
|                 |  |                  | reduction of nitrates to nitrogen  | colorless                         | pink-red                     |
|                 |  |                  |  | pink                              | Zn / 5 min<br>colorless      |
| TRP             | L - tryptophane                                | 0.2              | indole production (TRyptophane)  | JAMES / immediate                 |                              |
|                 |  |                  |  | colorless<br>pale green / yellow  | pink                         |
| GLU             | D-glucose                                      | 1.92             | fermentation (GLUcose)   | blue to green                     | yellow                       |
| ADH             | L-arginine                                     | 1.92             | Arginine DiHydrolase   | yellow                            | orange / pink / red          |
| URE             | urea   | 0.76             | UREase   | yellow                            | orange / pink / red          |
| ESC             | esculin<br>ferric citrate                      | 0.56<br>0.072    | hydrolysis ( $\beta$ - glucosidase) (ESCulin)                                  | yellow                            | grey / brown / black         |
| GEL             | gelatine<br>(bovine origin)                    | 0.6              | hydrolysis (protease) (GELatin)  | no pigment<br>diffusion           | diffusion of<br>black pigmen |
| PNPG            | 4-nitrophenyl- $\beta$ D-<br>galactopyranoside | 0.22             | $\beta$ - galactosidase (Para-NitroPhenyl -<br>$\beta$ D- Galactopyranosidase) | colorless                         | yellow                       |
| [GLU]           | D-glucose                                      | 1.56             | assimilation (GLUcose)   | transparent                       | opaque                       |
| [ARA]           | L-arabinose                                    | 1.4              | assimilation (ARABinose)   | transparent                       | opaque                       |
| [MNE]           | D-mannose                                      | 1.4              | assimilation (ManNosE)   | transparent                       | opaque                       |
| [MAN]           | D-mannitol                                     | 1.36             | assimilation (MANnitol)  | transparent                       | opaque                       |
| [NAG]           | N-acetyl-glucosamine                           | 1.28             | assimilation (N-Acetyl-Glucosamine)  | transparent                       | opaque                       |
| [MAL]           | D-maltose                                      | 1.4              | assimilation (MALtose)   | transparent                       | opaque                       |
| [GNT]           | potassium gluconate                            | 1.84             | assimilation (potassium GlucoNate)   | transparent                       | opaque                       |
| [CAP]           | capric acid                                    | 0.78             | assimilation (CAPric acid)   | transparent                       | opaque                       |
| [ADI]           | adipic acid                                    | 1.12             | assimilation (ADIpic acid)   | transparent                       | opaque                       |
| [MLT]           | malic acid                                     | 1.56             | assimilation (MaLaTe)  | transparent                       | opaque                       |
| [CIT]           | trisodium citrate                              | 2.28             | assimilation (trisodium CITrate)   | transparent                       | opaque                       |
| [PAC]           | phenylacetic acid                              | 0.8              | assimilation (PhenylACetic acid)   | transparent                       | opaque                       |
| OX              | (see oxidase test package insert)              | -                | cytochrome oxidase   | (see oxidase test package insert) |                              |

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