KOSSAN LATEX INDUSTRIES (M) SDN. BHD. (169832-K)

Rev. 004; Doc: SDS-003-PFN-B



Safety Data Sheet

Section 1: Identification

Product

Non-sterile, Ambidextrous, Single Use, Powder Free, Synthetic Rubber Nitrile Examination Gloves, Blue Coloured.

Identifier

PFN-B

Intended use

A patient examination glove is a disposable device intended for medical purposes that is worn on the examiner's hand or finger to prevent contamination between patient and examiner.

Manufacturer: Kossan Latex Industries (M) Sdn. Bhd

Address : Lot 6129, Jalan Haji Abdul Manan, Batu 5¼, Off Jalan Meru,

41050 Klang, Selangor Darul Ehsan, Malaysia.

Tel : +603-3392 3088 Fax : +603-3392 2328

Emergency Phone number

Tel : +603-3392 3088

Recommended use

For use in medical and healthcare environment.

Single use only.

Restrictions on use

Do not use for handling organic solvents.

Section 2: Hazard(s) Identification

1) Biological hazards

• Components used in glove manufacturing may cause allergic reaction in some user.

2) Hazards related to the re-use of the medical device

- Re-use of the non-sterile single use medical glove, leading to
 - i) loss of product lot traceability,
 - ii) risk of contamination and infection increases due to improper cleaning process,
 - iii) damages / weakening by the cleaning processes leading to increased risk of holes and tear during re-use.

Homepage: http://www.kossan.com.my

Rev. 004; Doc: SDS-003-PFN-B



KOSSAN LATEX INDUSTRIES (M) SDN. BHD. (169832-K)

Section 3: Composition/Information on Ingredients

Acrylonitrile Butadiene Rubber, Sulfur, Rubber Accelerator, Zinc Oxide, Potassium Hydroxide, Titanium Dioxide, Blue Pigment.

Section 4: First-Aid Measures

None, as glove is non-hazardous.

Section 5: Fire-Fighting Measures

Extinguishing media: Foam, Carbon Dioxide and Water

Flash point : N/A

Section 6: Accidental Release Measures

None.

Section 7: Handling and Storage

- The product should be stored in a cool dry place avoiding exposure to direct sunlight, high temperature and humidity.
- After the package is open, direct sunlight, fluorescent lighting and devices generating ozone should be avoided, and the product should be kept away from X-ray apparatus.
- Do not store these gloves with organic solvents as these solvents can degrade the gloves.
- Ventilation is not necessary under normal storage condition.

Section 8: Exposure Controls/Personal Protection

None.

Section 9: Physical and Chemical Properties

- 1) Feature
 These rubber gloves meet the ASTM D 6319 standard.
- 2) Form Soft solid form
- 3) Color Blue

Homepage: http://www.kossan.com.my



KOSSAN LATEX INDUSTRIES (M) SDN. BHD. (169832-K)

Rev. 004; Doc: SDS-003-PFN-B

Section 10: Stability and Reactivity

Stability : Stable Condition to avoid : None Incompatible Material : None

Section 11: Toxicological Information

Animal safety test for potential skin irritation and delayed skin sensitization showed that the glove is neither a skin irritant nor potential skin sensitizer.

Section 12: Ecological Information* (non-mandatory)

N/A

Section 13: Disposal Considerations* (non-mandatory)

The gloves are biodegradable product where no special decommissioning or disposal is required.

Section 14: Transport Information* (non-mandatory)

No special transportation requirements.

Section 15: Regulatory Information* (non-mandatory)

FDA Establishment Registration No. : 3005632797

• FDA 510(k) Notification No. : N/A

• FDA Device Listing No. : N/A

Section 16: Other Information

Date of preparation: 30th April 2015

Homepage: http://www.kossan.com.my

^{*}Note: Since other Agencies regulate this information, OSHA will not be enforcing Sections 12 through 15 (29 CFR 1910.1200(g)(2)).